


Product Datasheet

HDAC7 Antibody (orb1253010)

Description	HDAC7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	HDAC7A (phospho Ser155) antibody was raised against a peptide sequence around phosphorylation site of Serine 155 (T-V-S (p)-E-P) derived from Human HDAC7A.
Target	HDAC7
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	103 kDa
Uniprot ID	Q8WUI4
NCBI	NP_056216.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...