


Product Datasheet

GRIN2B Antibody (orb1253012)

Description	GRIN2B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	GRIN2B (phospho Ser1303) antibody was raised against a peptide sequence around phosphorylation site of Serine 1303 (Q-H-S (p)-Y-D) derived from Human GRIN2B.
Target	GRIN2B
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	170 kDa
Uniprot ID	Q13224
NCBI	NP_000825.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from M...