


Product Datasheet

FADD Antibody (orb1253013)

Description	FADD Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	FADD (phospho Ser191) antibody was raised against a peptide sequence around phosphorylation site of Serine191 (R-S-G (p)-A-M) derived from Human FADD.
Target	FADD
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	Q13158
NCBI	NP_003815.1
Dilution Range	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...



Immunohistochemical analysis of paraffin...