
Product Datasheet

NCF1 Antibody (orb1253022)

Description

NCF1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

p47 phox (phospho Ser345) antibody was raised against a peptide sequence around phosphorylation site of Serine 345 (P-Q-S (p)-P-G) derived from Human Nuclear p47 phox.

Target

NCF1

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW


45 kDa

Uniprot ID
[P14598](#)
NCBI
[NP_000256.4](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of paraffin...