


Product Datasheet

MDC1 Antibody (orb1253026)

Description	MDC1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	MDC1 (phospho Ser513) antibody was raised against a peptide sequence around phosphorylation site of Serine 513 (E-R-S (p)-Q-A) derived from Human MDC1.
Target	MDC1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	226 kDa
Uniprot ID	Q14676
NCBI	NP_055456.2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...