

---

## Product Datasheet

### GAB2 Antibody (orb1253029)

**Description**

GAB2 Antibody


**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

IHC, WB

**Applications**

Immunohistochemical analysis of paraffin...

**Immunogen**

Gab2 (phospho Ser623) antibody was raised against a peptide sequence around phosphorylation site of Serine 623 (P-S-S (p)-P-S) derived from Human Gab2 .

**Target**

GAB2

**Preservatives**

 Antibody supplied in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

Store antibody at -20°C for up to one year.

**Note**

For research use only

**Application notes**

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

80 kDa

**Uniprot ID**
[Q9UQC2](#)
**NCBI**
[NP\\_536739.1](#)
**Dilution Range**

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

**Expiration Date**

12 months from date of receipt.