
Product Datasheet

ESR1 Antibody (orb1253031)

Description

ESR1 Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

 Immunohistochemical
analysis of paraffin...

Immunogen

Estrogen Receptor-alpha (phospho Ser102) antibody was raised against a peptide sequence around phosphorylation site of Serine102 (L-N-S (p)-V-S) derived from Human Estrogen Receptor-alpha.

Target

ESR1

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

53 kDa

Uniprot ID
[P03372](#)
NCBI
[NP_000116.2](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.