


Product Datasheet

TCF3 Antibody (orb1253034)

Description	TCF3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	E2A (phospho Thr355) antibody was raised against a peptide sequence around phosphorylation site of threonine 355 (P-S-T (p)-P-V) derived from Human E2A .
Target	TCF3
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	P15923
NCBI	NP_003191.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...