


Product Datasheet

APLP2 Antibody (orb1253048)

Description	APLP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	APLP2 (phospho Tyr755) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 755 (P-T-Y (p)-K-Y) derived from Human APLP2.
Target	APLP2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86 kDa
Uniprot ID	Q06481
NCBI	NP_001633.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...