
Product Datasheet

PTK2B Antibody (orb1253056)

Description

PTK2B Antibody

Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications


Immunohistochemical analysis of paraffin...

Immunogen

PYK2 (phospho Tyr579) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 579 (E-D-Y (p)-Y-K) derived from Human PYK2.

Target

PTK2B

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

116 kDa

Uniprot ID
[Q14289](https://www.uniprot.org/uniprot/Q14289)
NCBI
[NP_004094.3](https://www.ncbi.nlm.nih.gov/nuccore/NP_004094.3)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.