


Product Datasheet

ITCH Antibody (orb1253057)

Description	ITCH Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	ITCH (phospho Tyr420) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 420 (F-I-Y (p)-G-N) derived from Human ITCH .
Target	ITCH
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102 kDa
Uniprot ID	Q96J02
NCBI	NP_001244066.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...