
Product Datasheet

FOS Antibody (orb1253067)

Description

FOS Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical analysis of paraffin...

Immunogen

FOS (phospho Thr232) antibody was raised against a peptide sequence around phosphorylation site of threonine 232 (V-A-T (p)-P-E) derived from Human FOS.

Target

FOS

Preservatives

 Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

48 kDa

Uniprot ID
P01100
NCBI
NP_005243.1
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.