


Product Datasheet

LATS1, LATS2 Antibody (orb1253092)

Description	LATS1, LATS2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	LATS1/2 (phospho Thr1079/1041) antibody was raised against a peptide sequence around phosphorylation site of threonine 1079/1041 (E-F-T (p)-F-R) derived from Human LATS1/2 .
Target	LATS1, LATS2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	126 kDa
Uniprot ID	O95835 , Q9NRM7
NCBI	NP_004681.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...