


Product Datasheet

IRAK1 Antibody (orb1253093)

Description	IRAK1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	IRAK1 (phospho Ser376) antibody was raised against a peptide sequence around phosphorylation site of Serine 376 (Q-S-S (p)-M-V) derived from Human IRAK1.
Target	IRAK1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76 kDa
Uniprot ID	P51617
NCBI	NP_001560.2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...