


Product Datasheet

EPHA3, EPHA4 Antibody (orb1253099)

Description	EPHA3, EPHA4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	EPHA3/4/5 (phospho Tyr779/833) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 779/833 (E-A-Y (p)-T-T)/ (Q-D-V (p)-I-K) derived from Human EPHA3/4/5.
Target	EPHA3, EPHA4
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110 kDa
Uniprot ID	P29320
NCBI	NP_004429.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...