


---

## Product Datasheet

### CSNK1A1, CSNK1A1L Antibody (orb1253100)

**Description**

CSNK1A1, CSNK1A1L Antibody


**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

IHC, WB

**Applications**

Immunohistochemical analysis of paraffin...

**Immunogen**

CK-1alpha (phospho Tyr294) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 294 (Y-D-Y (p)-T-F) derived from Human CK-1alpha.

**Target**

CSNK1A1, CSNK1A1L

**Preservatives**

 Antibody supplied in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

Store antibody at -20°C for up to one year.

**Note**

For research use only

**Application notes**

 Western Blot: 1:500~1:1000,  
Immunohistochemistry: 1:50~1:100

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

37 kDa

**Uniprot ID**
[Q8N752](#), [P48729](#)
**NCBI**
[NP\\_001883.4](#)
**Dilution Range**

 Western Blot: 1:500~1:1000,  
Immunohistochemistry: 1:50~1:100

**Expiration Date**

12 months from date of receipt.