


Product Datasheet

Alk Antibody (orb1253102)

Description	Alk Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	ALK (phospho Tyr1096) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 1096 (P-N-Y (p)-C-F) derived from Human ALK.
Target	Alk
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	176 kDa
Uniprot ID	Q9UM73
NCBI	NP_004295.2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...