
Product Datasheet

RB1 Antibody (orb1253108)

Description

RB1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Retinoblastoma (phospho Ser608) antibody was raised against a peptide sequence around phosphorylation site of Serine 608 (Y-L-S (p)-P-V) derived from Human Retinoblastoma.

Target

RB1

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW


110 kDa

Uniprot ID
[P06400](#)
NCBI
[NP_000312.2](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of paraffin...