


Product Datasheet

PDPK1 Antibody (orb1253109)

| | |
|----------------------------|---|
| Description | PDPK1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Immunogen | PDK1 (phospho Tyr9) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 9 (A-L-Y (p)-D-A) derived from Human PDK1. |
| Target | PDPK1 |
| Preservatives | Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1 mg/mL |
| Storage | Store antibody at -20°C for up to one year. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 63 kDa |
| Uniprot ID | O15530 |
| NCBI | NP_002604.1 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of paraffin...