


Product Datasheet

TP53 Antibody (orb1253110)

Description	TP53 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	p53 (phospho Thr387) antibody was raised against a peptide sequence around phosphorylation site of threonine 387 (F-K-T (p)-E-G) derived from Human p53.
Target	TP53
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P04637
NCBI	NP_000537.3
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...