


Product Datasheet

IGF1R Antibody (orb1253111)

Description	IGF1R Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	IGF1R (phospho Tyr1346) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 1346 (Q-P-Y (p)-A-H) derived from Human IGF1R .
Target	IGF1R
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	154 kDa
Uniprot ID	P08069
NCBI	NP_000866.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...