


Product Datasheet

HSF1 Antibody (orb1253112)

Description	HSF1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	HSF1 (phospho Thr142) antibody was raised against a peptide sequence around phosphorylation site of threonine 142 (L-L-T (p)-D-V) derived from Human HSF1 .
Target	HSF1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	Q00613
NCBI	NP_005517.1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...