
Product Datasheet

HNRNPD Antibody (orb1253127)

Description

HNRNPD Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC, WB

Applications
Immunogen

hnRPD (phospho Ser83) antibody was raised against a peptide sequence around phosphorylation site of Serine 83 (N-S-S (p)-P-R) derived from Human hnRPD

Target

HNRNPD

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100, Immunofluorescence: 1:100~1:200

Isotype

IgG

Clonality

Polyclonal

MW


38 kDa

Uniprot ID
[Q14103](#)
NCBI
[NP_112738.1](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100, Immunofluorescence: 1:100~1:200

Expiration Date

12 months from date of receipt.



Immunofluorescence staining of methanol-...