


Product Datasheet

EPHA2, EPHA3, EPHA4 Antibody (orb1253131)

Description	EPHA2, EPHA3, EPHA4 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	IF, WB	
Immunogen	EPHA 2/3/4 (phospho Tyr588/596) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 588/596 (K-T-Y (p)-V-R)/ (R-T-Y (p)-V-D) derived from Human EPHA2/3/4.	
Target	EPHA2, EPHA3, EPHA4	Immunofluorescence staining of methanol...
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	Store antibody at -20°C for up to one year.	
Note	For research use only	
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200	
Isotype	IgG	
Clonality	Polyclonal	
MW	130 kDa	
Uniprot ID	P29317 , P29320 , P54764	
NCBI	NP_004422.2	
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200	
Expiration Date	12 months from date of receipt.	