

Product Datasheet

CHEK2 Antibody (orb1253137)

Description	CHEK2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Chk2 (phospho Thr383) antibody was raised against a peptide sequence around phosphorylation site of threonine383 (M-R-T (p)-L-C) derived from Human Chk1.
Target	CHEK2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	O96017
NCBI	NP_009125.1
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of methanol-...