


Product Datasheet

BAX Antibody (orb1253140)

Description	BAX Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Bax (phospho Thr167) antibody was raised against a peptide sequence around phosphorylation site of threonine167 (P-G-T (p)-P-T) derived from Human Bax.
Target	BAX
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q07812
NCBI	NP_620116.1
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of methanol-...