
Product Datasheet

PEA15 Antibody (orb1253143)

Description

PEA15 Antibody

Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications


Immunohistochemical analysis of paraffin...

Immunogen

PEA-15 (phospho Ser104) antibody was raised against a peptide sequence around phosphorylation site of Serine 104 (I-P-S (p)-A-K) derived from Human PEA-15.

Target

PEA15

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

19 kDa

Uniprot ID
[Q15121](#)
NCBI
[NP_003759.1](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.