

Product Datasheet

CHEK1 Antibody (orb1253145)

Description

CHEK1 Antibody


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical analysis of paraffin...

Immunogen

Chk1 (phospho Ser296) antibody was raised against a peptide sequence around phosphorylation site of Serine 296 (I-F-S (p)-N-L) derived from Human Chk1.

Target

CHEK1

Preservatives

 Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

45 kDa

Uniprot ID
[O14757](#)
NCBI
[NP_001107593.1](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.