

Product Datasheet

ARRB1 Antibody (orb1253158)

Description

ARR1 Antibody


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical analysis of paraffin...

Immunogen

Arrestin 1 (phospho Ser412) antibody was raised against a peptide sequence around phosphorylation site of Serine 412 (T-G-S (p)-P-Q) derived from Human Arrestin 1.

Target

ARR1

Preservatives

 Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

50 kDa

Uniprot ID
[P49407](#)
NCBI
[NP_004032.2](#)
Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.