


Product Datasheet

NTRK1 Antibody (orb1253164)

Description	NTRK1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Trk A (phospho Tyr791) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 791 (P-V-Y (p)-L-D) derived from Human Trk A.
Target	NTRK1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140 kDa
Uniprot ID	P04629
NCBI	NP_002520.2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...