


Product Datasheet

Fbxl3 Antibody (orb1253276)

Description	Fbxl3 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	FBXL3 antibody was raised against a peptide sequence around aa.1~5 (M-K-R-G-G) derived from mouse FBXL3 .
Target	Fbxl3
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q8C4V4
NCBI	NP_056637.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...