

Product Datasheet

CuNVgp5 Antibody (orb1253306)

Description	CuNVgp5 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	RNA silencing suppressor antibody was raised against a peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.
Target	CuNVgp5
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P15184
NCBI	NP_040957.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.