



Product Datasheet

FAM129A Antibody (orb1253313)

| | |
|----------------------------|---|
| Description | FAM129A Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Immunogen | Niban antibody was raised against a peptide sequence around aa. 912~916 (K-E-G-E-G) derived from Human Niban. |
| Target | FAM129A |
| Preservatives | Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1 mg/mL |
| Storage | Store antibody at -20°C for up to one year. |
| Note | For research use only |
| Application notes | Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200 |
| Clonality | Polyclonal |
| MW | 150 kDa |
| Uniprot ID | Q9BZQ8 |
| NCBI | NP_443198.1 |
| Dilution Range | Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of lysed extracts ...



Immunofluorescence staining of methanol-...