



Product Datasheet

FAM129B Antibody (orb1253318)

Description	FAM129B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Niban-like protein (Ab-712) antibody was raised against a peptide sequence around aa.710~714 (Q-V-S-S-P) derived from Human Niban-like protein.
Target	FAM129B
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	Q96TA1
NCBI	NP_001030611.1
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...



Immunofluorescence staining of methanol-...