


Product Datasheet

PIK3R1 Antibody (orb1253435)

Description	PIK3R1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	PI3K P85 (phospho-Tyr467) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 467 (L-Y (p) -E-E-Y) derived from Human PI3K P85.
Target	PIK3R1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	P27986 , Q92569 , O00459
NCBI	NP_852556.2 , NP_005018.1 , NP_001107644.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...