



Product Datasheet

SHC1 Antibody (orb1253443)

Description	SHC1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Shc1 (Phospho-Tyr349) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 349 (H-Q-Y (p) -Y-N) derived from Human Shc1.
Target	SHC1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Clonality	Polyclonal
MW	46, 52 kDa
Uniprot ID	P29353
NCBI	NP_001123512.1
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...



Immunofluorescence staining of methanol-...