


Product Datasheet

Camk2a Antibody (orb1253448)

Description	Camk2a Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	CaMKII (Phospho-Thr286) antibody was raised against a peptide sequence around phosphorylation site of threonine 286 (Q-E-T (p) -V-D) derived from Human CaMKII.
Target	Camk2a
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	P11798
NCBI	NP_001273738.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...