


Product Datasheet

EGFR Antibody (orb1253455)

Description	EGFR Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	EGFR (phospho-Tyr1110) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 1110 (P-V-Y (p) -H-N) derived from Human EGFR.
Target	EGFR
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	175 kDa
Uniprot ID	P00533
NCBI	NP_005219.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...