


Product Datasheet

SP1 Antibody (orb1253459)

Description	SP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	SP1 (Phospho-Thr739) antibody was raised against a peptide sequence around phosphorylation site of threonine 739 (T-A-T (p) -P-S) derived from human SP1 .
Target	SP1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	P08047
NCBI	NP_612482.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...