


Product Datasheet

LCK Antibody (orb1253481)

Description	LCK Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Lck (Phospho-Tyr394) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 394 (N-E-Y (p) -T-A) derived from Human Lck.
Target	LCK
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P06239
NCBI	NP_001036236.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...