



Product Datasheet

TP53 Antibody (orb1253493)

Description	TP53 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	p53 (Phospho-Ser33) antibody was raised against a peptide sequence around phosphorylation site of serine 33 (V-L-S (p) -P-L) derived from Human p53.
Target	TP53
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	P04637
NCBI	NP_000537.3
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...



Immunofluorescence staining of methanol-...