


Product Datasheet

NBN Antibody (orb1253500)

Description	NBN Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	p95/NBS1 (Phospho-Ser343) antibody was raised against a peptide sequence around phosphorylation site of serine 343 (S-L-S (p) -Q-G) derived from Human p95/NBS1.
Target	NBN
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	95 kDa
Uniprot ID	O60934
NCBI	NP_002476.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...