


Product Datasheet

HNF4A Antibody (orb1253501)

Description	HNF4A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	HNF4a (Phospho-Ser304) antibody was raised against a peptide sequence around phosphorylation site of serine 304 (L-R-S (p) -Q-V) derived from Human HNF4 α .
Target	HNF4A
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:1000
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	P41235
NCBI	NP_000448.3
Dilution Range	Western Blot: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...