


Product Datasheet

JUN Antibody (orb1253502)

Description	JUN Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	c-Jun (Phospho-Tyr170) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 170 (P-V-Y (p) -A-N) derived from Human c-Jun.
Target	JUN
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P05412
NCBI	NP_002219.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...