


Product Datasheet

HDAC4 HDAC5 HDAC9 Antibody (orb1253696)

Description	HDAC4 HDAC5 HDAC9 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220) antibody was raised against a peptide sequence around phosphorylation site of serine 246/259/220 (T-A-S (p) -EP) derived from Human HDAC4/HDAC5/HDAC9.
Target	HDAC4 HDAC5 HDAC9
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Clonality	Polyclonal
MW	140, 124 kDa
Uniprot ID	P56524 , Q9UKV0 , Q9UQL6
NCBI	NP_055522.1 , NP_001015053.1 , NP_006028.2
Dilution Range	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...



Immunohistochemical analysis of paraffin...