



Product Datasheet

NFKBIA Antibody (orb1253760)

Description	NFKBIA Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	IκB-α (Phospho-Tyr42) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 42 (E-E-Y (p) -E-Q) derived from Human IκappaB-α.
Target	NFKBIA
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P25963
NCBI	NP_065390.1
Dilution Range	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extract from He...



Immunohistochemical analysis of paraffin...