



Product Datasheet

HDAC8 Antibody (orb1253780)

Description	HDAC8 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	HDAC8 (Phospho-Ser39) antibody was raised against a peptide sequence around phosphorylation site of serine 39 (R-A-S (p) -M-V) derived from Human HDAC8.
Target	HDAC8
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Immunohistochemistry: 1:50~1:100, Immunofluorescence: 1:100~1:200
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	Q9BY41
NCBI	NP_001159890.1
Dilution Range	Immunohistochemistry: 1:50~1:100, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunofluorescence staining of methanol-
...