
Product Datasheet

**rev, p19, Human immunodeficiency virus 1, HIV
Antibody (orb1253849)**

Description

rev, p19, Human immunodeficiency virus 1, HIV

Species/Host

Gallus

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

WB

Immunogen

1-116

Target

rev, p19, Human immunodeficiency virus 1, HIV

Preservatives

Phosphate-Buffered Saline. No preservatives added.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

HIV-Rev antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

HIV-Rev antibody can be used for the detection of HIV-Rev by Western Blot.

Clonality

Polyclonal

MW

13.1 kDa (calculated)

Uniprot ID
[P04618](#)
NCBI
[AAK08487.1](#)
Dilution Range

HIV-Rev antibody can be used for the detection of HIV-Rev by Western Blot.

Expiration Date

12 months from date of receipt.

