
Product Datasheet

PTPN13 Antibody (orb1254006)

Description

PTPN13 Antibody

Species/Host

Gallus

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

212-399

Target

PTPN13

Preservatives

Phosphate-Buffered Saline. No preservatives added.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PTPN13 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Protein-tyrosine phosphatase antibody can be used for the detection of Protein-tyrosine phosphatase by Western Blot.

Clonality

Polyclonal

MW


14.2 kDa (calculated)

Uniprot ID
[Q12923](#)
NCBI
[NP_006255.1](#)
Dilution Range

Protein-tyrosine phosphatase antibody can be used for the detection of Protein-tyrosine phosphatase by Western Blot.

Expiration Date

12 months from date of receipt.



E coli-derived fusion protein as test an...