
Product Datasheet

PPFIA1 Antibody (orb1254014)

Description

PPFIA1 Antibody

Species/Host

Gallus

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications


E coli-derived fusion protein as test an...

Immunogen

310-360

Target

PPFIA1

Preservatives

Phosphate-Buffered Saline. No preservatives added.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PPFIA1 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

PTPRF interacting protein α 1 isoform b variant antibody can be used for the detection of PTPRF interacting protein α 1 isoform b variant by Western Blot.

Clonality

Polyclonal

MW

77.8 kDa (calculated)

Uniprot ID
[Q13136](#)
NCBI
[NP_003617.1](#)
Dilution Range

PTPRF interacting protein α 1 isoform b variant antibody can be used for the detection of PTPRF interacting protein α 1 isoform b variant by Western Blot.

Expiration Date

12 months from date of receipt.